

Marquis™

A Fresh Scented Granular Presoak Concentrate

FEATURES

- Uniquely designed canister
- Specially formulated powder
- Alkaline / surfactant blend

BENEFITS

- Safe, convenient and simple to use
- Removes tough flatware soils
- Penetrates soils quickly

RECOMMENDED USES

For presoaking silver and stainlessware.

SPECIFICATIONS

Color : Blue

Scent/Odor : Citrus

Character : Granular powder

Foam : Low

pH use dilution : 9.5 - 11.0

Phosphorus content : 3.4%

Approvals : cRc



IMPORTANT

Read and understand the label and Safety Data Sheet (SDS). These provide a description of hazards, precautions, personal protective equipment, first aid, 24 Hour Emergency lines, storage data and other safe handling information. The SDS is provided with each container shipped. It is also available on our website or from the distributor.

PRODUCT DIRECTIONS

EQUIPMENT TO USE

Recommended warewash dispensing system.

DIRECTIONS TO USE

1) Dispense product solution into bus tub. Consult your sanitation specialist for proper use dilution. 2) Place stainless ware and/or silverware in solution for several minutes. 3) Remove, rinse and wash in dishmachine in normal fashion. **CANISTER SET UP:** 1) Remove seal from canister. 2) Place the canister into the Freedom System™ specialties concentrate dispenser, open end down. Use handle when transporting canister. 3) DO NOT remove concentrate canister until completely empty. 4) To remove empty canister, lift it from the dispenser carefully by the handle and turn upright to avoid spills. Dispose of canister immediately. DO NOT reuse empty canister. 5) Place new concentrate canister into dispenser. NOTE: Dispenser WILL NOT activate unless a canister is in unit.

PACKAGING

This product is available in 4/7 pound cases (5898964).

PROline™