EPA Reg. No. 1839-86-7546 EPA Est. No. 7546-WI-1 (WT); 875-KY-1 (KY)

Shurguard Plus



Quaternary Based Sanitizer for Institutional Use Technical Data Sheet

Shurguard Plus is designed for use in hospitals, nursing homes, hotels, food service establishments, restaurants and bars where disinfection, sanitization and deodorization is of prime importance. When used as directed, this product is formulated to disinfect hard, non-porous, inanimate environmental surfaces: floors, walls, glazed porcelain, glazed ceramic tile, plastic surfaces, bathrooms, cabinets, tables and chairs.

Specifications	
Color	Red
Scent/Odor	None added
Character	Clear, non-viscous liquid
Foam	Moderate
pH Use Dilution	7.0 - 8.0
Phosphorus Content	0.0%
Approvals	cRc, NSF D2
Registrations	EPA
Pack Sizes	
057602	4/1 gallon case
3296521	5 gallon pail



Directions For Use

Equipment: Jug or pail pump, wall pump or faucet feeder, recommended automatic dispensing system, trigger sprayer, mop, sponge, cleaning cloth, bowl mop, rubber gloves, eye protection and protective clothing.

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. This product is not for use on critical and semi-critical medical device surfaces.

DISINFECTION:

For visibly soiled areas, a pre-cleaning step is required. Prepare a fresh solution for each use. To disinfect inanimate, hard non-porous surfaces apply use-solution with mop, cloth, sponge, low pressure coarse sprayer or hand pump trigger spray so as to wet all surfaces thoroughly. Allow to remain visibly wet for 10 minutes, then remove excess liquid. Surfaces which may contact food must be rinsed thoroughly with potable water after use. For low pressure coarse spray and hand pump trigger sprayer applications, spray 6-8 inches from the surfaces; rub with brush, sponge

or cloth. Do not breathe spray.

NOTE: For spray applications, cover or remove all food particles.

Hospital Disinfection and Virucidal Activity - Add 3.5 fl. oz. of this product per 5 gal of water for hospital disinfection and virucidal activity.

SARS-CoV-2 Activity - At 2 fl. oz. per gal. this product is an effective virucide on hard, non-porous environmental surface against SARS-CoV-2 after 30 seconds.

FOOD CONTACT SURFACE SANITIZATION DIRECTIONS:

To sanitize pre-cleaned and potable water-rinsed, non-porous food contact surfaces, prepare a 200-400 ppm active quaternary solution by adding 1-2 fl. oz. of this product to 4 gal of water. immerse in a 200 to 400 ppm active quaternary solution for at least 60 seconds, making sure to immerse completely. Remove, drain the use-solution from the surfaces and air dry. Do not rinse.

Safety Reminder

Read and understand the label and Safety Data Sheet (SDS). These provide a description of hazards, precautions, personal protective equipment, first aid, 24 Hour Emergency lines, storage data and other safe handling information. The SDS is available on our website at www.uschemical.com.

Check with U S Chemical regarding product availability and registration status in your state. See product Reference Sheet or container label for more information.

